



- Insulated stainless steel cabinet
- Frost resistant down to (-40 °C)
- Additional frontal heating element (-20 °C to -40 °C)
- Fan for ventilation
- 2 x LED lights (*meaning defibrillator is visible in the dark*)
- Designed to be easy to open, even under conditions of heavy snow and frost.
- Future-proof! (*Software can be upgraded*)

FUNCTIONAL DESCRIPTION:

If the door of the cabinet is opened...	The SIREN immediately sounds. <i>(if the door of the cabinet is closed once again within 3 minutes the siren will stop automatically)</i>
Once the door is closed...	The cabinet is once again ready for operation.

TECHNICAL DESCRIPTION:

Door:	The door is fitted with thermopane glass and marked with international AED labels/symbols.
Light:	2 x LED lights in the cabinet, such that the defibrillator is visible in the dark.
Siren:	Siren in the cabinet. (100 dBA / 1m)
Fan:	Ventilation fan, starts automatically when temperature exceeds +40 °C.
Frost-proof cabinet:	Thermostatically-controlled heating element. Frost resistant down to -40 °C
Voltage:	Must be connected to a 230V (116Wh) electrical supply. Outdoor cable (3 x 1.5 mm ²) Cat II, IE IEC 61010-1
Overvoltage category:	2
Protection rating:	IP43.
Ambient temperature:	-40 / +40°C
Relative humidity:	0 – 95 %
Environmental pollution degree:	2
Altitude:	Max. 2000 masl. (<i>meters above sea level</i>)
Design norm:	IEC61010 – with deviations in insulation fire class and service switch.
Insulation – fire class:	E / F (EN13501)
External dimensions:	W x H x D = 44 x 46 x 27 cm (<i>incl. door handles</i>)
Internal dimensions:	W x H x D = 33 x 38 x 20 cm
Weight:	Gross weight incl. packaging 18.5 kg.

MAINTENANCE / CLEANING:

We recommend cleaning the outer surface of the cabinet at least 1-2 times per year.
Wipe using detergent and a damp cloth or brush and rinse with clean water. Finally wipe down the cabinet with a dry cloth.

If the surface of the cabinet becomes tarnished, this will be due to deposits which can easily be cleaned off.
For difficult deposits, use a scouring pad. **Never use wire wool.**