

microlife®

A partner for people. For life.



reddot award 2018
winner

3_{sec}

Automatic
measurement



NC₂₀₀

Non Contact Thermometer

- Touch-free measurement in only 3 seconds
- Gentle and hygienic
- Auto measurement with distance control
- 30 Data memories (stored with date and time)
- Fever alarm
- Multifunctional (body, ambient temperature)
- °C - °F switchable

Key features

Clinically
tested

Clinically tested

Several studies (conducted in Austria) prove the high accuracy and reliability of Microlife's infrared thermometers.



Auto measurement with distance control

CN Pat. ZL 201410095361.7; Others are patent pending

This thermometer takes automatically a measurement, when the device detects the distance is within 5 cm. Switch the device on and point to measure.



Guidance system for selfmeasurement

The self-indicator light on the back shows the user that the device is at the right distance.



User-friendly and hygienic

Fast and easy measurement without touching the person.

- Very easy handling
- Gentle monitoring, even on a sleeping child
- Blue tracking light for easy measurements at night

3in1




Multifunctional and precise

Accurately measures surface temperature of the following:

- Milk in a baby's bottle
- Surface of a baby's bath
- Ambient temperature


Measuring range and accuracy: 0 - 100.0 °C;
±0.2 °C between 35.0 ~ 42.2 °C



Silent Glow™ Technology – Fever alarm

The colour-illuminated fever alarm:

- Green display: temperatures up to 37.4 °C
- Red display: temperatures higher than 37.4 °C



Illuminated display with date and time

Large illuminated display with date and time.



Sound Control (On/Off)

Silent operation especially for measurements at night.