

microlife®

A partner for people. For life.



3_{sec}

with Bluetooth®
technology



NC150 BT

Non Contact Thermometer

- Touch-free measurement in only 3 seconds
- Gentle and hygienic
- Data transfer with Bluetooth® technology
- 30 Data memories (stored with date and time)
- Fever alarm
- Multifunctional (body, ambient temperature)
- °C – °F switchable

Key features

**Clinically
tested**

Clinically tested

Several studies (conducted in Austria) prove the high accuracy and reliability of Microlife's infrared thermometers.

3 in 1




Multifunctional and precise

Accurately measures surface temperature of:

- Milk in a baby's bottle
- Surface of a baby's bath
- Ambient temperature

Measuring range and accuracy: 0 – 100.0 °C;
± 0.2 °C between 35.0 ~ 42.2 °C



Silent Glow™ Technology – Fever alarm

The colour-illuminated fever alarm:

- Green display: temperatures up to 37.4 °C
- Red display: temperatures higher than 37.4 °C



User-friendly and hygienic

Fast measurement without touching the person.

- Gentle monitoring, even on a sleeping child
- Blue tracking light for easy measurements at night




Microlife Connected Health+ App

- Helps to check and analyse your health
- Using Bluetooth® Smart for easy data transfer
- Compatibility:
Android™ OS 5 or newer / iPhone iOS 9.0 or newer



Bluetooth® technology

Instant temperature monitoring on your smartphone.



Illuminated display with date and time

Large illuminated display with date and time.



Sound Control (On/Off)

Silent operation especially for measurements at night.